

KingBreathing Bags

anaesthesia breathing bags



Key Benefits

Comfortable

- Thin yet durable design
- Five standard sizes offer a bag to fit any patient
- The hourglass-shaped Contour breathing bag reduces hand fatigue

Gripable

- Tri-pleat design offers a better fit

Dependable

- Thin and durable surface maximizes clinician feel and ability to assess patient lung activity
- The hourglass-shaped Contour breathing bag efficiently delivers required gas volume

King Breathing Bags

King Breathing bags include a complete range of latex-free bags in sizes 0,5 ltr, 1,0 ltr, 2,0 ltr and 3,0 ltr with different types of connectors in compliance with ISO regulations. The Rainbow Breathing bags are a colourful option for

infants and children. In addition to pediatric anesthesia breathing bags, the Contour® breathing bags with an hourglass shape provide a better fit and grip.

KingBreathing Bags

anesthesia breathing bags

1/2 Liter Size

Description	Product #
OPKS* Racing Red	R2510
OPKS Sunshine Yellow	Y2510
OPKS Bouncin' Blue	B2510
OPKS Tickled Pink	K2510
OPKS Party Purple	P2510
OPKS Variety Pack	V2510
Green	20905

2 Liter Size

Description	Product #
OPKS Racing Red	R10902
OPKS Sunshine Yellow	Y10902
OPKS Bouncin' Blue	B10902
OPKS Tickled Pink	K10902
OPKS Party Purple	P10902
OPKS Variety Pack	V10902
Green	20902



Rainbow Breathing Bags®

*OPKS indicates Operation: Kid Stuff®. Rainbow breathing bags are a component of Operation: Kids Stuff products. All breathing bags are packaged 30 per case. All breathing bags are latex free.

1 Liter Size

Description	Product #
OPKS Racing Red	R10901
OPKS Sunshine Yellow	Y10901
OPKS Bouncin' Blue	B10901
OPKS Tickled Pink	K10901
OPKS Party Purple	P10901
OPKS Variety Pack	V10901
Green	20901

3 Liter Size

Description	Product #
Green	20903

Contour®

Description	Product #
Contour, Green, 2L	30902Z

Materials

Part	Material
Bag	Neoprene
Bushing	TPO (POE)
Tape	Vinyl Tape (PVC with DEHP)

Product weight

Size	Weight
1/2L	18.4g
1L	25.2g
2L	37.6g
3L	56.5g

Sales Company

Ambu A/S
 Baltorpbakken 13
 DK - 2750 Ballerup
 Denmark
 Tel. +45 7225 2000
 Fax +45 7225 2053
 www.ambu.com

Manufacturer

King Systems
 15011 Herriman Boulevard
 Noblesville, IN 46060
 Tel. 317 776 6823
 Tel. 800 642 5464
 Fax 317 776 6827
 www.Kingsystems.com



CE US: Rx only