

DATASHEET

Unique manikin helps to prepare students to respond to a more typical cardiac arrest event with a victim who is elderly and overweight.

Palpable anatomical landmarks including sternum, rib cage and substernal notch.

The manikin allows head tilt and chin lift to open the airways properly. The airway can be manipulated to simulate airway obstruction or choking situations.

Sanitary System

The individual mouth/nosepieces are easy to clean and the disposable airway system with a non-re-breathing valve prevents cross-contamination and is changed once a day.

Bariatric CPR Manikin



KEY BENEFITS

- More typical cardiac arrest victim for realistic CPR training
- Elderly and overweight appearance
- Airway opens with head tilt and chin lift
- Airway obstruction / choking simulation
- Sanitary system with mouth/nosepieces and disposable lung/airway system

ORDERING SPECIFICATIONS

Item no.

Bariatric CPR Manikin (Fat old Fred)

LF03750U	3 mouth/nosepieces, 3 disposable lung/airway systems, carry bag
LF03751U	Lung/Airway Systems Bariatric CPR Manikin (Fat old Fred), disposable (set of 24)
LF03752U	Mouth/Nosepieces Bariatric CPR Manikin (Fat old Fred), (set of 10)

Ambu

Ambu A/S
Baltorpbakken 13
2750 Ballerup
Denmark
T +45 72 25 20 00
F +45 72 25 20 50
ambu.com