

# SINGLE-USE LIKE YOU'VE NEVER SEEN IT BEFORE

Ambu® aScope™ 5 Broncho combines the excellent manoeuvrability and next-level imaging needed in the bronchoscopy suite with the sterility and efficiency of the single-use concept. The result? Consistently high performance, sterility, and the flexibility you need for a complex and demanding workflow. It is compatible with the full-HD resolution Ambu® aBox™ 2 displaying and processing unit, which can easily be connected to external monitors.




aScope 5 Broncho delivers consistently high performance and excellent manoeuvrability, supported by a new rotational function, and a smart functional design.

## EXCELLENT MANOEUVRABILITY

### Access difficult-to-reach areas of the bronchial tree:

The combination of the 195°/195° up and down bending angles, distal tip angulation range, and flexibility of the insertion tube helps you access difficult-to-reach areas in the bronchial tree. You will notice the excellent manoeuvrability, not least when working with an endotherapy instrument.

**120° left/right rotation** is designed to enable you to adjust the position of the distal end in order to ease navigation and insertion of endotherapy instruments.





Ambu® aScope™ 5 Broncho

## NEXT-LEVEL IMAGING

High-resolution camera chip and full-HD displaying and processing unit (12.8" screen): aScope 5 Broncho and aBox 2 work together to deliver superior imaging.

### Excellent image quality:

Advanced image processing, 120° field-of-view and 3-100 mm depth-of-field HD system with advanced adaptive auto control give you the bright and clear image needed for the procedure.

### Enhanced visibility:

Ambu's proprietary Advanced Red Contrast (ARC) enhancement algorithm is intended to improve the visibility of red colour tones in the image without altering the natural appearance.

### Easy to take, review and document images and videos

At the touch of a button on the bronchoscope, you can zoom in and out, capture images and videos and review them post-procedure with patients or colleagues. Compatibility with PACS servers and electronic medical record software makes documentation smooth and efficient.

## CONTACT US TO ARRANGE A DEMO

It's quick and simple for us to set up aScope 5 Broncho and aBox 2 at your place of work so you can get it in your hands and see for yourself.

**Ambu** FOREVER  
FORWARD

Ambu A/S  
Baltorpbakken 13  
DK - 2750 Ballerup  
Denmark  
Tel.: +45 7225 2000  
ambu.com

